
Originalbetriebsanleitung  
Translation of the original operating instruction

**STERIPUR 306 / KMA 305**  
**2/2 Wege-Membranventil**  
**2/2 way diaphragm valve**

---

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Translation of the original operating instructions

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STERIPUR 306 / KMA 305

### Table of contents

|        |                                                          |    |
|--------|----------------------------------------------------------|----|
| 1.     | Introduction.....                                        | 29 |
| 1.1.   | Information about the operating instructions.....        | 30 |
| 2.     | Safety.....                                              | 31 |
| 2.1.   | Safety Information.....                                  | 31 |
| 2.2.   | Hazard classification .....                              | 31 |
| 2.3.   | Intended use.....                                        | 32 |
| 2.4.   | Misuse .....                                             | 32 |
| 2.5.   | General safety provisions .....                          | 32 |
| 2.6.   | Residual risks .....                                     | 33 |
| 2.7.   | Responsibilities of the operator .....                   | 34 |
| 2.7.1. | Appointment and instruction of responsible persons ..... | 34 |
| 2.7.2. | Information obligation.....                              | 34 |
| 2.8.   | Description of target groups.....                        | 34 |
| 3.     | Transport / Storage / Disposal .....                     | 35 |
| 3.1.   | Delivery.....                                            | 35 |
| 3.1.1. | Scope of delivery.....                                   | 35 |
| 3.1.2. | Notification about damage .....                          | 35 |
| 3.2.   | Transport .....                                          | 35 |
| 3.3.   | Storage.....                                             | 36 |
| 3.4.   | Disposal .....                                           | 36 |
| 3.5.   | Return.....                                              | 36 |
| 4.     | Technical data .....                                     | 37 |
| 4.1.   | General technical data.....                              | 37 |
| 4.2.   | Operating conditions .....                               | 37 |
| 4.2.1. | Operating temperatures .....                             | 37 |
| 4.2.2. | Operating pressure .....                                 | 38 |
| 4.3.   | Safety related data .....                                | 39 |
| 4.4.   | Dimensions .....                                         | 40 |
| 4.5.   | Type plate .....                                         | 41 |
| 5.     | Installation.....                                        | 42 |
| 5.1.   | Tools required.....                                      | 42 |
| 5.2.   | Structure .....                                          | 43 |
| 5.3.   | Assembly & disassembly of the actuator .....             | 44 |
| 5.3.1. | Assembly.....                                            | 44 |
| 5.3.2. | Disassembly .....                                        | 45 |

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---

|        |                                                                |    |
|--------|----------------------------------------------------------------|----|
| 5.4    | Installation .....                                             | 46 |
| 5.4.1. | Self-draining of 2/2-way valves for aseptic applications ..... | 46 |
| 5.4.2. | Connection butt weld ends .....                                | 46 |
| 5.5.   | Implementation .....                                           | 47 |
| 5.6    | Operation of the integrated accessory .....                    | 48 |
| 5.6.1. | Optical position indicator .....                               | 48 |
| 5.6.2. | Integrated travel stop .....                                   | 48 |
| 5.6.3. | Integrated stroke limiter .....                                | 49 |
| 5.6.4. | Integrated fixing device of the handwheel .....                | 49 |
| 6.     | Maintenance .....                                              | 50 |
| 6.1.   | Wear parts .....                                               | 50 |
| 6.2.   | Diaphragm change .....                                         | 52 |
| 6.2.1. | Change of one-piece diaphragms .....                           | 52 |
| 6.3.   | Cleaning .....                                                 | 53 |
| 7.     | Accessories .....                                              | 53 |
| 8.     | Troubleshooting .....                                          | 54 |
| 9.     | EC Declaration of conformity .....                             | 55 |

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EC declaration of conformity acc. Directive 2014/68 / EU

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Bad Rappenau, 09.04.2024

Uwe Rutsch  
(Head of Development)

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 1. Introduction

In case you have any queries about the device, please contact our customer services stating the serial number:

#### **SED Flow Control GmbH**

Am Schafbaum 2

D-74906 Bad Rappenau

PO-Box 1306

D-74900 Bad Rappenau

Phone: +49(0)7264/921-0

Fax: +49(0)7264/921-21

E-Mail: [info-sed-de@samsongroup.com](mailto:info-sed-de@samsongroup.com)

Web: [sed.samsongroup.com](http://sed.samsongroup.com)

WEE Registration:

SED Flow Control GmbH

Registration ID: 30963364

These operating instructions describe the technical status of the device at delivery.

It is valid for standard versions. For special versions please contact our customer services.

Any subsequent changes at the device are not considered in these operating instructions.

© 2025

The present operating instructions document is protected by copyright. Without written approval by the manufacturer it may not be translated, reproduced (mechanically or electronically) or disclosed to third parties, either in full or in parts.

The German version of these operating instructions is deciding in case of doubt, misunderstandings, or translation mistakes!

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 1.1. Information about the operating instructions

#### Safe operation

The operating instructions contain important information for safe and correct installation of the device. Compliance with that helps preventing hazards, avoiding repair costs and downtime, and increasing reliability and operating life of the device.

#### Target groups

The contents of these operating instructions are intended for installation and maintenance staff.

#### Reading the operating instructions

All persons who work at the device must read the operating instructions in order to be familiar with correct handling and operation. The operating instructions contain important safety information! Non-compliance with such information may cause hazardous situations. The operating instructions must be read and understood.

#### Copyright

The present operating instructions document is protected by copyright. Without written approval by SED Flow Control GmbH it must not be reproduced, distributed, or unauthorizedly used for competitive purposes, either in full or in parts.

All rights for drawings and other documents and any power of disposal are solely with SED Flow Control GmbH; this also applies for cases of application for intellectual property rights.

#### Declaration of conformity

The device complies with the basic requirements of the applicable European directives. Conformity was proved.

#### Technical changes

We constantly strive for further development and improvement of our devices. The data contained in this edition correspond to the state of the art at the time of printing. All technical details with regard to information and figures within the operating instructions are subject to change.

#### Warranty

The present printed document does not contain any guarantee commitments. For that purpose, please see our general terms and conditions of delivery and payment. The prerequisite for statutory warranty is the intended use of the device in compliance with the specified operating conditions.

The warranty only applies to absence of defectiveness of the device and its components. No liability will be assumed for consequential damage occurring through failure or malfunction of the device, irrespective of their nature.

#### Information on the internet

Instructions and data sheets for the device are provided on the internet at:

[sed.samsongroup.com](http://sed.samsongroup.com)

Translation of the original operating instructions  
STERIPUR 306 / KMA 305

## 2. Safety

### 2.1. Safety Information

#### Warning

The operating instructions contain important safety information!  
Non-compliance with such information may cause hazardous situations.  
The operating instructions must be read and understood.

### 2.2. Hazard classification

#### Danger

Indicates imminent danger!  
Non-compliance will result in death or severe injuries.

#### Warning

Indicates a potentially hazardous situation!  
Non-compliance may result in severe injuries.

#### Attention

Indicates a possible hazard!  
Non-compliance may result in moderate or slight injuries.

#### Notice

Warning from property damage!  
Non-compliance may cause damage to the device or the facility.

Indicates important additional information, tips and recommendations

Indicates important reference to information in these operating instructions and in other documentation.

Indicates an operating step that must be performed.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 2.3. Intended use

The diaphragm valves STERIPUR 306 and KMA 305 are designed for the installation in pipelines and are used to control and regulate media.

Only expose the devices to media to which the body and seals are resistant. Intended use also includes compliance with the installation and maintenance conditions specified by the manufacturer as well as the general rules of technology. Any other use is considered improper. The manufacturer is not liable for any damage resulting from this; the user bears the risk.

#### **Warning**

The relevant and generally accepted safety rules apply when planning the operation, as well as operating the device. Appropriate measures must be taken to prevent unintentional actuation or undue impairment. Please note that planners, plant builders or operators are generally responsible for the positioning and installation of the electrical position indicators.

### 2.4. Misuse

- Do not supply the system with flammable or aggressive medias.
- Do not apply mechanical loads on the housing (e.g. positioning of objects on the device, using it as a step).
- Do not carry out any external modifications on the housings of devices. Do not paint any parts of the housing.
- Use the device only with components or spare parts, which are approved from SED Flow Control GmbH.
- The standard version of the diaphragm valve must not be used in potentially explosive areas. For this purpose, please use the suitable version and use it only in approved regions!

### 2.5. General safety provisions

- The operating instructions by the manufacturer are binding for installation, programming, maintenance and transport.
- It is the operator's responsibility to complement the safety provisions by special instructions specific for local operating conditions, if required.
- The operating instructions and all safety-relevant rules and regulations must be kept in a safe place.
- The operating instructions and safety information must be available in full and in easily legible condition.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### Accident prevention rules

In addition to the operating instructions and the binding rules and regulations for accident prevention applicable in the country of use and the place of operation, the approved technical rules for safe and professional working must be observed.

### Before beginning with your work

Before beginning with your work, please obtain information about first-aid and rescue options (emergency, doctor, fire- and rescue services). Please obtain information about locations and operation of fire extinguishers and about fire alarm and firefighting options on site. Retain the unit against inadvertent operation.

### During work

During operation, wear suitable work clothes.  
Do not carry out activities that may impair safety.  
Only operate the device in safe and functional condition.

### Provisions for protection of the environment

For any work at or with the device, the provisions for environmental protection must be complied with.

## 2.6. Residual risks

### Danger

Danger of injury through high pressure!

Pressure acting on lines and valves may cause severe injuries!

Before disconnecting any lines and valves, the pressure must be switched off and the lines must be vented.

At no time overstep the maximum operation or control pressure!

Danger through electric voltage!

Electrical voltage may cause severe injuries or death!

Before working at the facility, voltage must be switched off and secured against unauthorized activation!

Attend to valid accident prevention and safety guidelines for electrical devices!

### Warning

Protect the facility from inadvertent operation!

After shutdown the facility a controlled restart must be guaranteed!

### Attention

Danger through hot device surface while continuous operation!

Hot device surface may cause burning - or danger of fire!

Keep the device away from easy burning materials and do not touch with bare hands!

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 2.7. Responsibilities of the operator

Devices that are not in unobjectionable condition may cause personal injury and property damage.

It is the operator's responsibility to only operate the device in unobjectionable condition.

Hazard zones that occur between the device and any equipment on the customer's site must be secured by the operator. The operator must ensure sufficient illumination during work at the device.

#### 2.7.1. Appointment and instruction of responsible persons

- Only employ staff who has received safety instructions.
- Responsibilities and competencies of staff with regard to installation, programming and maintenance must be clearly defined.
- Regularly monitor the staff for safe and responsible working methods and compliance with the operating instructions

#### 2.7.2. Information obligation

It is the obligation of the operator of the device to make these operating instructions available to all persons who work with the device at any time.

Before use of the device, all persons must have read and fully understood the operating instructions.

### 2.8. Description of target groups

The contents of these operating instructions are intended for various target groups. The level of knowledge/information the respective target group must have is defined below.

All target groups must have read and fully understood the contents of these operating instructions.

Installation, programming and maintenance staff must

- Be 18 years of age.
- Have profound education and professional training.
- Be trained with regard to conduct in case of a fault.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 3. Transport / Storage / Disposal

#### 3.1. Delivery

##### 3.1.1. Scope of delivery

- Immediately after receipt of the delivery, please verify that the contents are not damaged, and the type and scope of the delivery complies with the data stated on the delivery note.
- By checking the order number, please verify that the goods received comply with the versions you ordered.
- In case of irregularities, please contact us immediately.

##### 3.1.2. Notification about damage

- Damage resulting from inadequate packaging or during transport must be communicated to the shipping company, the insurance, and the supplier immediately after delivery of the consignment.

#### 3.2. Transport

##### Notice

If the guidelines are not reserved, the device may be damaged!

- The device must be transported in impact-resistant packaging.
- During transport, the device must be protected against humidity and dirt.
- The approved ambient temperature of  $-10\text{ °C}$  to  $+55\text{ °C}$  must not be exceeded through the impact of heat or cold.

## Translation of the original operating instructions

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#### Danger

Danger of injury after reassembly!

Please check, if there are any damages and ensure, that the mounting is correct, especially loosened erection screws!

In order to keep a device that is not used over a longer period functional, a few items must be considered:

- Please store the device in its original packaging. The storage location must be dry and clean.
- Storage temperature must be between  $-10\text{ °C}$  and  $+50\text{ °C}$
- Lock up the control air ports with protecting caps.
- Secure the device against unauthorized operation.
- Please do not store in abrasive surrounding area.

### 3.4. Disposal

#### Protect the environment!

The individual components and the packaging must be properly disposed of according to the types of materials. Please note the current rules for waste disposal!

#### Notice

Please check all media contacting parts for contamination, if necessary dispose them according to the current rules of waste disposal and environmental regulations!

### 3.5. Return

#### Notice

Due to legal regulations for the protection of the environment and personnel, it is necessary that the return declaration is completely filled out and signed and enclosed with the shipping documents. Only if this declaration is completed in full will the return be processed!

- Clean the unit.
- Request a return declaration from SED.
- Return only with completely filled out return declaration.
- Otherwise
  - no credit note or no
  - completion of the repair
  - but a disposal with costs

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 4. Technical data

#### 4.1. General technical data

|                                    |                                                                                                                                                                                                                                                                                                            |                                                         |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| Series                             | STERIPUR 306                                                                                                                                                                                                                                                                                               | KMA 305                                                 |
| Size                               | DN 8 – DN 20                                                                                                                                                                                                                                                                                               |                                                         |
| Diaphragm size MA                  | MA 10                                                                                                                                                                                                                                                                                                      |                                                         |
| Actuator suitable for valve bodies | Two-Way bodies<br>Welded valve configuration<br>T-bodies<br>Multiport bodies<br>Tank bottom bodies                                                                                                                                                                                                         |                                                         |
| End connection valve body          | Butt weld ends acc. DIN EN 10357 (formerly DIN 11850 row 1 & 2)<br>Butt weld ends acc. DIN 11866 row B (formerly ISO 1127)<br>Butt weld ends acc. DIN 11866 row C (formerly ASME BPE / ASTM 269)<br>SMS 3008, JIS G 3447<br><br>Clamps, flanges and threads acc. Customer specifications or customer order |                                                         |
| Material actuator                  | Handwheel: Stainless steel 1.4408<br>Bonnet: Stainless steel 1.4408                                                                                                                                                                                                                                        | Handwheel: PARA GF 50<br>Bonnet: Stainless steel 1.4408 |
| Material valve body                | 1.4435 / 316L Investment cast<br>1.4435 / 316L forged<br>1.4435 / 316L forged Fe < 0,5% (Basler Norm)<br>Other materials acc. Customer specifications or customer order                                                                                                                                    |                                                         |
| Material diaphragm                 | EPDM, PTFE-TM / EPDM<br>Other materials acc. Customer specifications or customer order                                                                                                                                                                                                                     |                                                         |
| Installation position              | Any, actuator preferably upwards (except tank valves)                                                                                                                                                                                                                                                      |                                                         |

### 4.2. Operating conditions

#### 4.2.1. Operating temperatures

#### Warning

**Do not exceed or deceed the permitted temperatures!**

The permissible temperatures depend on the materials used. Please note that only the lowest value is permissible for the maximum medium temperature! Conversely, the highest value is permissible for the minimum medium temperature. At increased temperatures, the maximum permissible operating pressure may be reduced!

#### Permissible media temperature for actuators

160°C dependent on application.

STERIPUR 306: Autoclavable.

KMA 305: Suitable for a limited number of cycles of autoclaving

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### Permissible temperature for diaphragms

| Diaphragm material | Code | Media temperature |       | Sterilization temperature |
|--------------------|------|-------------------|-------|---------------------------|
|                    |      | Min.              | Max.  |                           |
| EPDM               | 28   | -10 °C            | 90 °C | 150 °C, 60 min            |
| PTFE-TM / EPDM     | 51   | -10 °C            | 90 °C | 150 °C, 60 min            |

### Permissible temperature for valve bodies

| Material  | Media temperature |        |
|-----------|-------------------|--------|
|           | Min.              | Max.   |
| St. Steel | -10 °C            | 160 °C |

## 4.2.2. Operating pressure

### Warning

Danger of bursting at excess pressure! Do not exceed the permitted operating pressure!  
At increased temperatures, the maximum permissible operating pressure may be reduced!

### Maximal permissible pressure, unidirectional, static pressure

|                |                   |
|----------------|-------------------|
| Diaphragm      | MA10<br>DN 8 - 20 |
| EPDM           | 10 bar            |
| PTFE-TM / EPDM | 10 bar            |

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 4.3. Safety related data

The diaphragm valves STERIPUR 306 / KMA 305 are in accordance to

- 3-A Standards Section 54-02
- TA-Luft / VDI 2440 / VDI 3479
- Pressure Equipment Directive 2014/68/EU (declaration of conformity, see attachment)

An ATEX version according to Directive 2014/34/EU is available for use in potentially explosive areas.

#### 4.4. Dimensions


Figure 1: Mounting dimensions


## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### Butt weld tube end standards

| Stutzenstandard                                           |       |     |        |     | ISO 1127<br>DIN 11866<br>Reihe B | DIN 11850 |                                 | DIN<br>Vor- zugs-<br>reihe | ASTM 269<br>ASME BPE DIN<br>11866<br>Reihe C | BS O.D.<br>4825                                                        | SMS 3008  | JIS G 3447            | JIS G 3459      |
|-----------------------------------------------------------|-------|-----|--------|-----|----------------------------------|-----------|---------------------------------|----------------------------|----------------------------------------------|------------------------------------------------------------------------|-----------|-----------------------|-----------------|
|                                                           |       |     |        |     |                                  | Reihe 1   | Reihe 2<br>DIN 11866<br>Reihe A |                            |                                              |                                                                        |           |                       |                 |
| Code                                                      |       |     |        |     | 40                               | 41        | 42                              | 39                         | 45 <sup>1</sup>                              | 94                                                                     | 49        | 97                    | 98 <sup>3</sup> |
| DN                                                        | NPS   | MA  | L(min) | L1  | ød x s                           | ød x s    | ød x s                          | ød x s                     | ød x s                                       | ød x s                                                                 | ød x s    | ød x s                | ød x s          |
| Ventiltyp handgesteuert 205 / 206 / 290 / 297             |       |     |        |     |                                  |           |                                 |                            |                                              | Ventiltyp fremdgesteuert 190 / 207 / 217 / 327                         |           |                       |                 |
| 4                                                         | -     | 8   | 20     | 72  | -                                | -         | -                               | 6x1,0                      | -                                            | -                                                                      | -         | -                     | -               |
| 6                                                         | -     | 8   | 20     | 72  | -                                | -         | 8x1,0 <sup>2</sup>              | 8x1,0                      | -                                            | -                                                                      | -         | -                     | 10,5x1,2        |
| 8                                                         | 1/4   | 8   | 20     | 72  | 13,5x1,6                         | -         | 10x1,0 <sup>2</sup>             | 10x1,0                     | 6,35x0,89                                    | -                                                                      | -         | -                     | 13,8x1,65       |
| 10                                                        | 3/8   | 8   | 20     | 72  | -                                | 12x1,0    | 13x1,5                          | -                          | 9,53x0,89                                    | -                                                                      | -         | -                     | -               |
| 15                                                        | 1/2   | 8   | 20     | 72  | -                                | -         | -                               | -                          | 12,7x1,65                                    | 12,7x1,2                                                               | -         | -                     | -               |
| Ventiltyp handgesteuert 289 / 295 / 397                   |       |     |        |     |                                  |           |                                 |                            |                                              | Ventiltyp fremdgesteuert 188 / 195 / 317 / 327 / 392 / 394             |           |                       |                 |
| 8                                                         | -     | 10  | 25     | 108 | 13,5x1,6                         | -         | -                               | -                          | -                                            | -                                                                      | -         | -                     | -               |
| 10                                                        | 3/8   | 10  | 25     | 108 | 17,2x1,6                         | 12x1,0    | 13x1,5                          | -                          | 9,53x0,89 <sup>3</sup>                       | -                                                                      | -         | -                     | 17,3x1,65       |
| 15                                                        | 1/2   | 10  | 25     | 108 | 21,3x1,6                         | 18x1,0    | 19x1,5                          | 18x1,5                     | 12,7x1,65                                    | 12,7x1,2                                                               | -         | -                     | 21,7x2,1        |
| 20                                                        | 3/4   | 10  | 25     | 108 | -                                | -         | 23x1,5                          | 22x1,5                     | 19,05x1,65                                   | 19,05x1,2                                                              | -         | -                     | -               |
| Ventiltyp handgesteuert 905 / 907 / 982 / 985 / 995 / 997 |       |     |        |     |                                  |           |                                 |                            |                                              | Ventiltyp fremdgesteuert 327 / 385 / 395 / 402 / 407 / 417 / 495 / 592 |           |                       |                 |
| 15                                                        | -     | 25  | 25     | 120 | 21,3x1,6                         | 18x1,0    | 19x1,5                          | -                          | 12,7x1,65 <sup>3</sup>                       | -                                                                      | -         | -                     | 21,7x2,1        |
| 20                                                        | 3/4   | 25  | 25     | 120 | 26,9x1,6                         | 22x1,0    | 23x1,5                          | -                          | 19,05x1,65                                   | -                                                                      | -         | -                     | 27,2x2,1        |
| 25                                                        | 1     | 25  | 25     | 120 | 33,7x2,0                         | 28x1,0    | 29x1,5                          | 28x1,5                     | 25,4x1,65                                    | -                                                                      | 25,0x1,2  | 25,4x1,2              | -               |
| 32                                                        | 1 1/4 | 40  | 25     | 153 | 42,4x2,0                         | 34x1,0    | 35x1,5                          | -                          | 31,75x1,65 <sup>3</sup>                      | -                                                                      | 33,7x1,2  | 31,8x1,2 <sup>3</sup> | -               |
| 40                                                        | 1 1/2 | 40  | 25     | 153 | 48,3x2,0                         | 40x1,0    | 41x1,5                          | -                          | 38,1x1,65                                    | -                                                                      | 38,0x1,2  | 38,1x1,2              | -               |
| 50                                                        | 2     | 50  | 30     | 173 | 60,3x2,0                         | 52x1,0    | 53x1,5                          | -                          | 50,8x1,65                                    | -                                                                      | 51,0x1,2  | 50,8x1,5              | -               |
| 65                                                        | 2 1/2 | 50  | 30     | 173 | -                                | -         | -                               | -                          | 63,5x1,65                                    | -                                                                      | 63,5x1,6  | 63,5x2,0 <sup>3</sup> | -               |
| 65                                                        | 2 1/2 | 80  | 25     | 216 | 76,1x2,0                         | -         | 70x2,0                          | -                          | 63,5x1,65                                    | -                                                                      | 63,5x1,6  | 63,5x2,0 <sup>3</sup> | -               |
| 80                                                        | 3     | 80  | 30     | 254 | 88,9x2,3                         | -         | 85x2,0                          | -                          | 76,2x1,65                                    | -                                                                      | 76,1x1,6  | 76,3x2,0              | -               |
| 100                                                       | 4     | 100 | 30     | 305 | 114,3x2,3                        | -         | 104x2,0                         | -                          | 101,6x2,11                                   | -                                                                      | 101,6x2,0 | 101,6x2,1             | -               |
| 125                                                       | 5     | 150 | 30     | 450 | 139,7x2,6                        | -         | 129x2,0                         | -                          | -                                            | -                                                                      | -         | -                     | -               |
| 150                                                       | 6     | 150 | 30     | 450 | 168,3x2,6                        | -         | 154x2,0                         | -                          | 152,4x2,77                                   | -                                                                      | -         | -                     | -               |

### 4.5. Type plate

The nameplate is attached to the device in the form of an adhesive label and contains important information about operation.

You can determine a detailed breakdown of the match code online at [sed.samsongroup.com](http://sed.samsongroup.com) using our product configurator.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 5. Installation

#### **Danger**

Danger of injury through high pressure!

Pressure acting on lines and valves may cause severe injuries!

Before disconnecting any lines and valves, the pressure must be switched off and the lines must be vented.

Danger of injury through electric shock!

Shut off current and secure against unauthorized activation.

#### **Warning**

Danger of injury in case of incorrect assembly!

Incorrect assembly may cause severe injuries.

Assembly and installation must only be performed by authorized expert staff.

Use suitable tools only!

Secure plant against unauthorized activation!

Guarantee controlled restart after shut down!

#### 5.1. Tools required

The tools required for mounting and installation of the device are not included in the scope of delivery:

Only use functional, suitable and safe tools!

5.2. Structure


Figure 2: Exploded view of complete valve

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### Functional description

The manually operated 2/2-way diaphragm valve STERIPUR 306 consists of a stainless steel bonnet and a stainless steel handwheel. KMA 305 consists of a stainless steel bonnet and a plastic handwheel. The design is characterized by a rising handwheel with visual indicator and stroke indication. It also has an integrated seal adjuster. KMA 305 has additionally a lockable handwheel. An integrated stroke limiter is available for both valves as an option.

### 5.3. Assembly & disassembly of the actuator

#### Danger

Danger of injury!

Check the device for correct assembly before implementation! Also check the tightening torque!

Only use functional, suitable and safe tools! Assembly and disassembly must only be performed by authorized expert staff.

#### 5.3.1. Assembly

- Check the technical data of the valve and the materials and ensure suitability for the respective application. The valve must be designed for the operating conditions of the pipeline system, e.g. medium, medium concentration, temperature and pressure, as well as the respective ambient conditions.
- Shut down the plant or plant section and secure it against being switched on again.
- Depressurize the system or system section.
- Completely drain the plant or plant section.
- Allow the plant or plant section to cool down until the evaporation temperature of the medium has fallen below to avoid scalding.
- Decontaminate, rinse and ventilate the system or system section properly.
- Move actuator to open position.
- Position the actuator (pos. 1) correctly on the valve body (pos. 2), see Fig. 2. Make sure that the alignment of the diaphragm web and the valve body web match!
- Position the fastening screws and tighten them hand-tight.
- Close pneumatic actuator to 80%.
- Tighten the screws (pos. 5) or nuts crosswise with the appropriate torque, see table below.

Table 1: Tightening torques screws - stainless steel body

| Membran-<br>größe<br><br><i>diaphragm<br/>size</i> | Gewinde-<br>größe<br><br><i>thread<br/>size</i> | Anzahl<br><br><i>No. of<br/>bolts</i> | Edelstahl                   |            |
|----------------------------------------------------|-------------------------------------------------|---------------------------------------|-----------------------------|------------|
|                                                    |                                                 |                                       | <i>stainless steel body</i> |            |
|                                                    |                                                 |                                       | PTFE                        | EPDM       |
|                                                    |                                                 |                                       | all                         | all        |
|                                                    |                                                 |                                       | Nm                          | Nm         |
| <b>MA10</b>                                        | <b>M5</b>                                       | <b>4</b>                              | <b>3-4</b>                  | <b>1-2</b> |

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

→ Check the valve for function and tightness

Take care for simultaneous compression of the diaphragm – about 10%, depending on diaphragm material and version. Especially PTFE-diaphragms have a setting behavior. Check the valve after the first thermal load latest and re-tighten the bolts if necessary.

### 5.3.2. Disassembly

#### **Danger**

**Danger of injury through high pressure!**

Pressure acting on lines and valves may cause severe injuries!

Before disconnection any lines and valves, the pressure must be switched off and the lines must be vented.

Secure plant against unauthorized activation!

**Danger of injury through media leakage (steam, hot medias, acid, lye)!**

**Important:** After disassembly, clean all parts from dirt (do not damage parts). Check parts for damage, replace if necessary (use only original SED parts).

- Consider suitable protective equipment according to the regulations of the plant operator.
- Shut down the plant or part of the plant.
- Secure against being switched on again.
- Depressurize the plant or plant section.
- Move valve to open position.
- Loosen fastening screws (pos. 5) crosswise.
- Remove actuator (pos. 1) incl. diaphragm (pos. 3) from valve body (pos. 2).
- Clean parts in case of contamination. Replace defective parts only with original spare parts from SED Flow Control.

For diaphragm assembly see chapter 6.2.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 5.4 Installation


Before installation, ensure that the diaphragm valve is suitable for the operating conditions. These include function, materials, operating pressure and working media in connection with operating temperature. Depending on the operating conditions and weight of the valve a suitable fixture has to be applied.

Take also care for aligning pipes.

Clean pipes before installation.

#### 5.4.1. Self-draining of 2/2-way valves for aseptic applications

The 2/2-way diaphragm valve can be installed independent of position and flow direction. For horizontal installation in self-draining position, the self-draining angle has to be followed:

The individual self-draining angle is marked on each valve body. The current revision of the self-draining angle table can be found in our catalog, which is available for download on our homepage: [www.sed.samsongroup.com](http://www.sed.samsongroup.com).

#### 5.4.2. Connection butt weld ends

Follow the laws, directives and regulations for welding which are valid in your country.

- Disassemble the valve before welding. See chapter 5.3.2.
- Weld in valve body and cool down. If necessary, weld in in self-draining position, see marking on valve body. Polish surface if needed.
- Assemble actuator and diaphragm on valve body, see chapter 5.3.1.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 5.5. Implementation

Ensure correct installation before implementation. Check tightness in clamping area of actuator, diaphragm and valve body. Also check the valve for function.

#### Warning

Aggressive chemicals! Chemical burns! Check the tightness of the media connections before start-up! Leak test only with suitable protective equipment.

#### Attention

Prevent against leakage! Provide protective measures against exceeding the maximum permissible pressure due to possible pressure surges (water hammer).

Before cleaning or commissioning the system:

For new equipment and after repairs, flush piping system with valve fully open (to remove harmful foreign matter).

#### Cleaning:

The operator of the plant is responsible for selecting the cleaning medium and carrying out the procedure.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 5.6 Operation of the integrated accessory

The diaphragm valve has an integrated stroke limiter with visual position indicator as standard. This makes it possible to set the partial stroke for fine dosing.


#### 5.6.1. Optical position indicator

The optical position indicator shows the actual valve position. Thus, it can be recognized whether the valve is open- or closed.

Valve position „close“

Valve position „open“

#### 5.6.2. Integrated travel stop

The integrated travel stop enables limiting closing position. This could be issued to guarantee a minimum flow or to protect the diaphragm from overloading when closing.

##### Adjust travel stop:

- Take the valve into the desired position. Please note that the maximum opening is up to 25%!
- Loosen the screw that secures the anti-rotation device.
- Lift up the anti-rotation device, providing traction between the spindle and the handwheel.
- Turn the handwheel down as far as it will go.
- Press the anti-rotation device into the gear and secure it with the screw.
- The travel stop has been adjusted.

An instruction video is also available under <https://youtu.be/PQbHyQjuZ9I?feature=shared>

##### Revoke travel stop:

- Bring valve in closed position.
- Loosen the screw that keeps the anti-rotation device fastened.
- Lift up the anti-rotation device, providing traction between the spindle and the handwheel.
- Turn the handwheel 720° in open direction.
- Press the anti-rotation device into the gear and secure it with the screw.
- Turn the handwheel in close direction until a resistance can be felt.
- Loosen the screw that keeps the anti-rotation device fastened.
- Lift up the anti-rotation device again.
- Turn the handwheel clockwise to the lower stop and then turn it back by 90°.
- Press the anti-rotation device into the gear and secure the screw with 2.5 Nm.
- The valve is now set back to the factory setting.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 5.6.3. Integrated stroke limiter

The integrated stroke limiter enables to adjust the opening stroke and to set the flow rate to a maximum value. To adjust the stroke limiter please follow the guidelines:

- Take the valve into the desired position. E.g. valve opening 75%!
- Loosen the screw that keeps the anti-rotation device fastened.
- Lift up the Anti-rotation device.
- By turning the anti-rotation device, move the stroke limiter nut to the lower end position.
- Press the Anti-rotation device into the gear and secure it with the screw with 2.5 Nm.
- The stroke limiter has been adjusted.

An instruction video is also available under <https://youtu.be/PQbHyQjuZ9I?feature=shared>


### 5.6.4. Integrated fixing device of the handwheel

The integrated fixing device of the hand wheel enables to fix the valve position. This function is only available for KMA 305.

To manipulate the fixing device please follow the guidelines:

- Take the valve into the desired position.
- The locking can only be operated when the mark on the handwheel matches the mark on the optical position indicator.

Lock the integrated fixing device:

16 different positions can be applied to the hand wheel during a full rotation.

An optional locking system to protect the fixing from unauthorized adjustment can be provided.

An instruction video is also available under <https://youtu.be/PQbHyQjuZ9I?feature=shared>

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 6. Maintenance

Maintenance is only allowed to be carried out by trained staff.

Diaphragm valves are low-maintenance parts. They have to be checked and maintained periodically. These periods are depending on working load and have to be defined according to valid rules.

#### Warning

Pressurized fittings! Risk of serious injury or death!

Only work on a depressurized system.

#### Attention

Hot system parts! Burns! Only work on cooled down plant.

#### Attention

Maintenance and servicing activities only by trained specialist personnel. SED accepts no liability for damage caused by improper handling or external influence. In case of doubt, contact SED before commissioning.

1. consider suitable protective equipment according to the regulations of the plant operator.
2. Shut down the plant or part of the plant.
3. secure against restarting.
4. depressurize the plant or plant section.

The operator must carry out regular visual inspections of the ball valves in accordance with the operating conditions and the hazard potential to prevent leaks and damage. The ball valve must also be disassembled and checked for wear at appropriate intervals.

#### 6.1. Wear parts

The diaphragm is a wearing part and is subject to natural wear.

The actuator of the diaphragm valve is low-maintenance, provided that the regulations of these operating instructions for use are observed.

The actuator, the medium diaphragm and the valve body are available as spare parts.

To obtain the relevant spare parts, quote the match code and serial number on the nameplate at [info-sed-de@samsongroup.com](mailto:info-sed-de@samsongroup.com) or contact your local representative at [sed.samsongroup.com](http://sed.samsongroup.com).

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### **Danger**

**Danger of injury!**

**Bonnet is under spring pressure!**

**Do not open the actuator! If actuator is damaged, replace the complete actuator!**

The lifetime of the diaphragm is dependent on following factors:

- Diaphragm material
- Media
- Operating temperature
- Operating pressure
- Control pressure
- Correct assembly and use

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

### 6.2. Diaphragm change

The diaphragm has to be controlled for wear after max. 100.000 ( $\geq$ MA80 50.000) switching operations. Depending on working load, the control period has to be reduced according to valid regulations and rules.

Particles in the working media can be abrasive – reduce control period accordingly.

Only use diaphragms which are suitable for your application.

Don't use higher control pressure than necessary. This helps to save the diaphragm.

#### Notice

Media leakage by wrong assembly!  
Ensure correct assembly before startup!

#### Danger

Danger of injury through high pressure!  
Pressure acting on lines and valves may cause severe injuries!  
Before disconnecting any lines and valves, the pressure must be switched off and the lines must be vented.

Protect the facility from inadvertent operation!  
Danger through media leakage (steam, hot media, acid, lye)

#### 6.2.1. Change of one-piece diaphragms

- Disassemble valve, see chapter 5.3.2.
- Bring valve in close position.
- Screw out diaphragm anticlockwise ( $\geq$ MA10).
- Insert compressor and position it correctly.
- Screw in diaphragm by hand up to the stop or respectively plug it in. Then position diaphragm by turning it back anticlockwise in a way that the wholes of the diaphragm and valve body fit together. Attention! Do not turn back more than 180°! Also take care for correct position of the diaphragm flap. If the diaphragm isn't assembled correctly, point loading can occur! This can lead to leakage and damage of the diaphragm.
- Bring valve in open position.
- Place actuator on the valve body. Take care for correct orientation of diaphragm and weir.
- Position screws and fasten hand-tight.
- Close manual actuator to 80%.
- Fasten screws or nuts crosswise.
- Tightening torques see chapter 5.3.1.
- Check the valve at the pass and to the outside for tightness and function. Take care for equal compression of the diaphragm – about 10% depending on material and version. Especially PTFE diaphragms show a setting behavior. Check the valve after the first thermal load latest and re-fasten screws if necessary.

## Translation of the original operating instructions

STERIPUR 306 / KMA 305

---

### 6.3. Cleaning

For external cleaning, commercially available cleaning agents can be used. Check the chemical resistance of all materials used in the valve, which could come in contact with the cleaning agent, before applying

For sterilization processes, check if the valve is suitable for the applied process.

### 7. Accessories

Fitting and functional accessories are available via [sed.samsongroup.com](https://sed.samsongroup.com) or contact us.

Translation of the original operating instructions  
STERIPUR 306 / KMA 305

## 8. Troubleshooting

| Fault                              | Potential cause                                                   | Solution                                                                               |
|------------------------------------|-------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Valve doesn't close or close fully | Media pressure too high                                           | Do not exceed maximum control pressure, see chapter 4.2.2.                             |
|                                    | Foreign matters between weir and diaphragm                        | Disassemble valve, remove foreign matters and check parts for damage. See Chapter 5.3. |
|                                    | Defect diaphragm                                                  | Change diaphragm, see chapter 6.2.ff                                                   |
|                                    | Integrated travel stop is still set                               | Adjust the travel stop acc. chapter 5.6 ff                                             |
| Valve doesn't open or open fully   | Wrong assembled diaphragm                                         | Disassemble actuator and assemble diaphragm correctly. See chapter 6.2.ff              |
|                                    | Stroke limiter is still set                                       | Adjust the stroke limiter acc. chapter 5.6 ff                                          |
| Media leakage at diaphragm flange  | Wrong assembled diaphragm                                         | Disassemble actuator and assemble diaphragm correctly. See chapter 6.2.ff              |
|                                    | Tightening torque of screw connection actuator -diaphragm too low | Check connection of actuator and valve body, see chapter 5.3.1.                        |
|                                    | Defect diaphragm                                                  | Change diaphragm, see chapter 6.2.ff                                                   |
|                                    | Defect actuator                                                   | Change actuator, see chapter 5.3.ff                                                    |
|                                    | Defect valve body                                                 | Change valve body                                                                      |
|                                    | Media pressure too high                                           | Do not exceed maximum control pressure, see chapter 4.2.2.                             |
|                                    | Operating temperature too high                                    | Do not exceed maximum operating pressure, see chapter 4.2.1.                           |
| Media leakage at pipe connection   | Defect sealing / sealing material                                 | Change sealing / sealing material                                                      |
|                                    | Improper installation                                             | Install valve body proper                                                              |
|                                    | Threaded ends loose                                               | Re-tighten connection                                                                  |
| Media leakage at valve body        | Defect valve body                                                 | Change valve body                                                                      |

# Translation of the original operating instructions

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A complete technical documentation in the sense of the directive is available. The operating instructions belonging to the product are available in the original version. The undersigned manufacturer hereby certifies that the design, manufacture and testing of these valves comply with the requirements of the Pressure Equipment Directive 2014/68/EU of the European Parliament and of the Council of 13 May 2014. The Pressure Equipment Directive is implemented in accordance with the AD 2000 regulations.

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Valves DN ≤ 25 comply with the Pressure Equipment Directive 2014/68/EU Art. 4 (3). Therefore, they must not be marked with a CE mark or the number of a notified body.

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## MARKETS

- Chemicals and petrochemicals
- Power and energy
- District heating and cooling, building automation
- General industry
- Industrial gases
- Food and beverages
- Metallurgy and mining
- Oil and gas
- Pharmaceuticals and biotechnology
- Marine equipment
- Water and wastewater
- Pulp and paper

## PRODUCTS

- Valves
- Self-operated regulators
- Actuators
- Valve accessories
- Signal converters
- Controllers and automation systems
- Sensors and thermostats
- Digital solutions

## SALES SITES

- More than 50 subsidiaries in over 40 countries
- More than 200 representatives

## PRODUCTION SITES

- SAMSON Germany, Frankfurt, established 1916  
Total plot and production area: 150,000 m<sup>2</sup>
- SAMSON France, Lyon, established 1962  
Total plot and production area: 23,400 m<sup>2</sup>
- SAMSON Turkey, Istanbul established 1984  
Total plot and production area: 11,053 m<sup>2</sup>
- SAMSON USA, Baytown, TX, established 1992  
Total plot and production area: 9,200 m<sup>2</sup>
- SAMSON China, Beijing, established 1998  
Total plot and production area: 10,138 m<sup>2</sup>
- SAMSON India, Pune district, established 1999  
Total plot and production area: 18,000 m<sup>2</sup>
- SAMSON AIR TORQUE, Bergamo, Italy  
Total plot and production area: 27,684 m<sup>2</sup>
- SAMSON CERA SYSTEM, Hermsdorf, Germany  
Total plot and production area: 14,700 m<sup>2</sup>
- SAMSON KT-ELEKTRONIK, Berlin, Germany  
Total plot and production area: 1,060 m<sup>2</sup>
- SAMSON LEUSCH, Neuss, Germany  
Total plot and production area: 18,400 m<sup>2</sup>
- SAMSON PFEIFFER, Kempen, Germany  
Total plot and production area: 35,400 m<sup>2</sup>
- SAMSON RINGO, Zaragoza, Spain  
Total plot and production area: 18,270 m<sup>2</sup>
- SAMSON SED, Bad Rappenau, Germany  
Total plot and production area: 10,370 m<sup>2</sup>
- SAMSON STARLINE, Bergamo, Italy  
Total plot and production area: 26,409 m<sup>2</sup>
- SAMSON VDH PRODUCTS, Netherlands
- SAMSON VETEC, Speyer, Germany  
Total plot and production area: 27,090 m<sup>2</sup>

SED Flow Control GmbH  
Am Schafbaum 2 · 74906 Bad Rappenau  
Telefon: +49 7264 921 0 · Telefax: +49 7264 921 21  
E-Mail: [info-sed-de@samsongroup.com](mailto:info-sed-de@samsongroup.com) · Internet: [sed.samsongroup.com](http://sed.samsongroup.com)